

qEvaGreen® 5X qPCR Mix with ROX, no ROX, ROX vial, capillary, or HRM

for research and *in vitro* use only

Recommended Reaction Mix

COMPONENT	CONC.	AMOUNT	FINAL CONC.
qEvaGreen®	5X	4 µL	1X
Primer, F ¹	10 pmol/µL	0.16 - 0.5 µL	80 - 250 nM
Primer, R ¹	10 pmol/µL	0.16 - 0.5 µL	80 - 250 nM
Template DNA ²		2 - 400 ng	0.1 - 20 ng/µL
Nuclease-free water		variable	
Total		20 µL	

¹ To maximize assay sensitivity use lowest concentration possible without compromising reaction efficiency. Doubling the reverse primer concentration may improve performance. Further optimization can be tested using final concentrations of 100 - 400 nM.
² Genomic DNA, plasmid DNA, or cDNA can be used as template. For cDNA up to 10% of the final reaction may be cDNA (i.e. for a 20 µL qPCR reaction, use up to 2.0 µL of undiluted cDNA).

Recommended Cycling Settings

CYCLE STEP	TEMP (°C)	TIME	CYCLES
Initial denaturation	95	15 min	1
Denaturation	95	15 sec	
Annealing	T _M - 4 ¹	20 sec	40
Extension	72	20 sec	

Recommended Light Cycler 480 Settings

	TEMP (°C)	ACQUISITION MODE	HOLD (hh:mm:ss)	RAMP RATE (°C/s)	ACQUISITION RATE (per °C)
				96-well	384-well
Initial Denaturation	95	None	00:15:00	4.4	4.8
Amplification	95	None	00:00:10	4.4	4.8
	T _M - 4 ¹	None	00:00:15	2.2	2.5
	72	Single	E ²	4.4	4.8
High Resolution Melt ³	95	None	00:01:00	4.4	4.8
	40	None	00:01:00	2.2	2.5
Cooling	65	Continuous	00:00:01	1	1
	95	Continuous	-	-	-
Cooling	40	Non	00:00:10	2.2	2.5
	20	-	-	-	-

¹ Set annealing temperature to be 4°C lower than T_M of primers or set it to 60-65°C initially.
² Set elongation hold time E = Amplicon Length/15 (i.e. use 30 sec for a 450 bp amplicon).
³ Use HRM cycle only for qARTA HRM mixes, else skip HRM settings.

Estimating primer melting temperature:

For primers of less than 25 nucleotides, T_M = 4 (G + C) + 2 (A + T), where G, C, A, T represent the number of respective nucleotides in the primer. For longer primers, use specialized software to calculate T_M.

Description

qEvaGreen® qPCR mix is a ready-to-use cocktail optimized containing all components necessary for amplification and detection of DNA in qPCR. Simply add nuclease-free water, template and primers.

- **HRM option.** HRM with ROX or without ROX mixes are optimized for High Resolution Melt.

EvaGreen® dye

EvaGreen® Dye is a green fluorescent nucleic acid dye replacing SYBR Green I, with added benefits:

- **Highly sensitive.** Exhibits lower PCR inhibition.

which allows higher dye concentrations that in turn produces greater qPCR signals.

- **Nonmutagenic and noncytotoxic.** Impermeable to cell membranes.

Select the standard settings for SYBR Green or FAM to use qEvaGreen mixes.

Eco-friendly shipping & storage

- Shipped at ambient temperature without ice, foam, or other wasteful packaging.
- Use within 2 weeks of arrival or store at -20°C.

qEvaGreen® Order Information

Cat. No.	Desc.	Vol.	No. Rxns 20 µL each	Mix Composition				Nuclease free H ₂ O
				qHot Taq	Final qEvaGreen® 5X Buffer	dNTPs dNTP, dCTP, & dTTP	BSA	
CGWR-01	1 mL	250					4 mL	
CGWR-02	with ROX	2 mL	500				8 mL	PRE-MIXED
CGWR-05	5 mL	1250					-	-
CGNR-01	no ROX	1 mL	250				4 mL	
CGNR-02	with ROX	2 mL	500	✓	WITH EVA GREEN® DYE	1.5 mM	8 mL	No
CGNR-05	5 mL	1250					-	-
CGRV-01	with ROX	1 mL	250				4 mL	
CGRV-02	with ROX	2 mL	500				8 mL	
CGRV-05	vial	5 mL	1250				-	-
CGC-02	capillary	2 mL	500				8 mL	
CGC-05	5 mL	1250					-	-
CGHWR-02	HRM with ROX	2 mL	500	✓	WITH EVA GREEN®	1.5 mM	8 mL	PRE-MIXED
CGHWR-05	with ROX	5 mL	1250				-	-
CGHNR-02	HRM no ROX	2 mL	500	✓	HRM BUFFER WITH EVA GREEN®	2.5 mM	8 mL	PRE-MIXED
CGHNR-05	with ROX	5 mL	1250				-	-

* If volume is designated, then that volume of 0.1 mM ROX is in a separate vial.

Some applications of this product are covered by patents issued or pending other than qARTa Bio, and may require a license which is not provided by the purchase of this product. User should obtain a patent license if appropriate.

Revision date: 2011.09.01

© 2011 qARTaBio, Inc.

T_ +1-310-326-1542

E_technical.support@qartabio.com

E_sales@qartabio.com

W_www.qartabio.com



RTA
B I O