

qTaq 5X PCR Master Mix with MgCl₂ and Ready-to-Load option

For *in vitro* use only

Protocol

1. Thaw and briefly centrifuge all components
2. Prepare your reactions using the Recommended Reaction Mix
3. Briefly centrifuge your samples
4. Run PCR using Recommended Cycling Settings

Recommended Reaction Mix		
COMPONENT	CONC.	AMOUNT
qTaq MM	5X	4 µL
Primer, F ¹	10 pmol/µL	0.2 - 0.6 µL
Primer, R ¹	10 pmol/µL	0.2 - 0.6 µL
Template DNA ²		2 - 400 ng
Nuclease-free water		variable
Total		20 µL

1 To maximize assay sensitivity use lowest concentration possible without compromising reaction efficiency. Doubling the reverse primer concentration may improve performance. Further optimization can be tested using final concentrations of 100 - 400 nM.

2 Further optimizations can vary depending on your DNA template.

Recommended Cycling Settings		
CYCLE STEP	TEMP (°C)	TIME
Initial denaturation	95	3 - 5 min
Denaturation	95	20 - 40 sec
Annealing	T _M - 4 ¹	30 - 60 sec
Extension	72	1 min/kb ²
Final extension	72	5 - 10 min

1 Set annealing temperature to be 4°C lower than T_M of primers.
2 Use amplicon length to optimize extension time: approx. 1 min per 1000 bases.

Estimating primer melting temperature:
For primers containing less than 25 nucleotides, T_M = 4 (G + C) + 2 (A + T), where G, C, A, T represent the number of respective nucleotides in the primer. If primers contain more than 25 nucleotides specialized software is recommended to calculate T_M.

qTaq 5X Master Mix Order Information

Cat. No.	Vol. (mL)	No. Pkts/20 µL each	Mix Composition						
			qTaq	5X Buffer	Final MgCl ₂ Conc.	Final dNTP Conc.	Blue Loading Dye	Yellow Loading Dye	Nuclease free H ₂ O
QTMML15-S	0.2 mL	50	✓	400 mM Tris-HCl	1.5 mM	2.00 µM each of dATP, dCTP, dGTP, dTTP	No	No	1 mL
QTMML15-01	1 mL	250		100 mM (NH ₄) ₂ SO ₄	2.5 mM	0.1% w/v Tween-20	Migration equivalent to 3.5 - 4.5 kb DNA fragment	Migration rate in excess of primers in 1% agarose gel.	4 mL
QTMMLR15-01	1 mL	250		100 mM (NH ₄) ₂ SO ₄	2.5 mM	0.1% w/v Tween-20	Migration rate in excess of primers in 1% agarose gel.	4 mL	
QTMMLR25-S	0.2 mL	50		100 mM (NH ₄) ₂ SO ₄	2.5 mM	0.1% w/v Tween-20	Migration rate in excess of primers in 1% agarose gel.	4 mL	
QTMMLR25-01	1 mL	250		100 mM (NH ₄) ₂ SO ₄	2.5 mM	0.1% w/v Tween-20	Migration rate in excess of primers in 1% agarose gel.	4 mL	

Description

qTaq 5X Master Mix is a ready-to-load cocktail for PCR. It contains all the components necessary for DNA amplification in PCR. Simply add nuclease-free water (supplied with the kit), template and primers.

Fidelity

qTaq fidelity is ~2.5 X 10⁶ errors per nucleotide incorporation event or ~4.0 x 10⁶ nucleotides incorporated before any error occurs.

Eco-friendly shipping & storage

- Shipped at ambient temperature without ice to reduce packaging waste
- Use within 2 weeks of arrival or store at -20°C

Some applications of this product are covered by patents issued to parties other than qARTa Bio, and may require a license which is not provided by the purchase of this product. User should obtain a patent license if appropriate.